

Cable Sets for Fixed Cable Installation



Installation sets for fixed cable installation, including C-Boxes for wall mounting in treatment and control rooms

Features

- ▶ Help to protect the sensitive connection cables of therapy dosimetry equipment against mechanical damage
- ▶ Let you avoid handling of long cables during routine work
- ▶ Help to protect the sensitive connectors of dosimetry cables

The cable installation sets are used for fixed installation of dosimetric connection cables between treatment rooms and control rooms. Two C-Boxes mounted to the walls in the treatment room and in the control room serve as cable terminal boxes. The dosimetry equipment is simply connected to the C-Boxes by means of short-length cables that are removed after use.

Two versions of cable sets for firm installation are available:

- ▶ Cable set for the connection of the TBA control unit installed in the control room,
- ▶ Cable set for the connection of the TBA control unit installed in the treatment room.

If other dosimetry equipment, like the MULTIDOS has to be connected, additional cable components are required, depending on the various applications. Please ask for our suggestions to achieve the optimum solution.

Ordering Information

Installation sets for control unit set up in control rooms:

L981062 TBA cable installation set including cable for field and reference detector, cable length 20 m

L981063 TBA cable installation set including cable for field detector, cable length 20 m

T16006.1.001 C-Box for wall mounting (2 units required)

Installation set for control unit set up in treatment rooms:

T25023 RS232 cable installation set