


# CTDI Set for CT Dosimetry



## *Accessories for determination of CTDI in computed tomography*

### Features

- ▶ For measurements acc. IEC 61223-3-5 and 60601-2-44
- ▶ CT ion chamber, sensitive length 10 cm or 30 cm
- ▶ Body phantom, 32 cm Ø, and head phantom, 16 cm Ø, both acrylic cylinders of 15 cm height
- ▶ Combined head and body phantom is available

The CT chamber T30009 is inserted for CTDI measurements into one of the five holes of the head or body phantom. Holes not used are filled by acrylic dummy plugs. Etched crosshairs on the phantoms allow exact alignment.

The CT chamber T30017 is used for dose length product (DLP) measurements free in air.

### Ordering Information

*CT chambers, connecting system BNT, TNC or M:*

- 30009 CT chamber, length 10 cm
- 30017 CT chamber, length 30 cm
- T40016 CT body measuring phantom
- T40017 CT head measuring phantom
- T40027 CT head and body phantom

### Options

- T16018 DIADOS CT adapter
- T40016.1.010 Carrying case for body phantom
- T40017.1.010 Carrying case for head phantom, DIADOS and CT chamber
- T40027.1.050 Carrying case for head and body phantom, DIADOS E and CT chamber

- ▶ UNIDOS E *page 90*
- ▶ DIADOS, DIADOS E *page 88*